

## Water Retention Centrifuge



Centrifuge with 4 receptacles for dewatering containers. Dewatering containers with screen and sump for centrifuged water.

## Features

- Centrifuge with 4 dewatering containers

## Specifications

- Centrifuge with 4 receptacles for dewatering containers
- 100 to 14,500 rpm
- dewatering container volume: 22 or 43 ml
- Weight: approx. 150g /container
- Optionally available with Voith system attachments (4 receptacles each with 4 small dewatering containers)
- In this version 16 small samples can be centrifuged simultaneously

## Test description

Wet pulp samples are put into the dewatering containers which are inserted into the receptacles in the centrifuge. After the programmed time of dewatering each sample is taken out and weighed to calculate the water-retention value.

## Applicable standards

ISO 23714  
SCAN C 62:00

## Technical data

Electricity: 230 V, 50/60 Hz  
Water supply: None  
Compressed air: None  
Dimensions (WxDxH): 350x400x300 cm  
Net weight: 34 kg  
Gross weight: 40 kg