

| N<br>O. | Equipment<br>description            | Rotation<br>Flask 20L  | Rotation Flask<br>30L               | Rotation Flask<br>50L                | Rotation Flask<br>100L               | Rotation Flask<br>200L               |
|---------|-------------------------------------|--|-------------------------------------|--------------------------------------|--------------------------------------|--------------------------------------|
|         |                                     | RE1-20G-<br>1HP  | RE1-30G-1HP                         | RE1-50G-1HP                          | RE1-100G-1HP                         | RE1-200G-1HP                         |
|         |                                     | Ascending  |                                     |                                      |                                      |                                      |
| 1       | Bath Lift                           | Hand Lift  | Hand Lift                           | Hydraulic Lift                       | Hydraulic Lift                       | Hydraulic Lift                       |
| 2       | Circulation<br>Pump (for<br>shower) | YES  | YES                                 | YES                                  | YES                                  | YES                                  |
| 3       | Dimension                           | 1110W*650<br>D*2150H   | 1515W*650D*2<br>200H                | 1940W*900D*2<br>900H                 | 2050W*900D*2<br>970H                 | 2355W*900D*3<br>520H                 |
| 4       | Condensate<br>Area                  | 0.6m <sup>2</sup>  | 1m <sup>2</sup>                     | 1.25m <sup>2</sup>                   | 2.3m <sup>2</sup>                    | 4.0m <sup>2</sup>                    |
| 5       | Rotation<br>Speed                   | 0 ~120 rpm   | 0 ~120 rpm                          | 12 ~120 rpm                          | 12 ~120 rpm                          | 0 ~ 60 rpm                           |
| 6       | Receiving<br>Flask                  | 10L  | 20L                                 | 20L                                  | 30L                                  | 50L                                  |
| 7       | RPM<br>indicator                    | Digital Display  |                                     |                                      |                                      |                                      |
| 8       | Vacuum<br>(Ultimate<br>Pressure)    | 133Pa(1torr)   | 133Pa(1torr)                        | 133Pa(1torr)                         | 133Pa(1torr)                         | 667Pa(5torr)                         |
| 9       | Temp.Rang<br>e Heating<br>Bath      | Digital Display: 20 ~ 90°C (Water Bath)/Digital: 20 ~ 200°C (Oil Bath) |                                     |                                      |                                      |                                      |
| 10      | Heater                              | Total 4kW<br>(main 2kW+<br>sub 2kW)                                    | Total 4kW<br>(main 2kW+<br>sub 2kW) | Total 10kW<br>(main 5kW+<br>sub 5kW) | Total 12kW<br>(main 6kW+<br>sub 6kW) | Total 12kW<br>(main 6kW+<br>sub 6kW) |
| 11      | Temp.<br>Control                    | PID from main heater, ON/OFF Control                                   |                                     |                                      |                                      |                                      |
| 12      | Overheat<br>Cut-off<br>Protection   | YES  |                                     |                                      |                                      |                                      |
| 13      | Power<br>Source                     | AC200V 3<br>phase,20A  | AC200V 3<br>phase,20A               | AC200V 3<br>phase,30A                | AC200V 3<br>phase,40A                | AC200V 3<br>phase,50A                |