

OE Series

Ultimate quality control



Our range of optical emission spectrometers for QA/QC and process control

HITACHI

The background of the page features several thick, curved, metallic-looking bands that sweep across the frame from the top right towards the bottom left. The bands are rendered in shades of grey and white, with highlights and shadows that give them a three-dimensional, industrial appearance. The overall composition is clean and modern, emphasizing precision and technology.

For fast, comprehensive metals quality analysis in a single affordable instrument, the OE series from Hitachi High-Tech delivers everything you need.

High performance and efficiency

The OE720 and OE750 are groundbreaking OES metals analyzers. Covering the complete spectrum of elements in metal, they also have low detection limits.

Tightening of industry regulations, complex supply chains and increased use of scrap as a base material means it's crucial for foundries and metals manufacturers to control tramp and trace elements in the lowest ppm range. Historically, OES analysis at this level was out of reach for many businesses. That's now changed with the OE series from Hitachi High-Tech.

These spark spectrometers allow you to analyze all main alloying elements and identify exceptionally low levels of tramp, trace and treatment elements in

metals, such as nitrogen in steel. Fast measurement times, high reliability and low operating costs mean the OE720 and OE750 are invaluable for everyday analysis and total quality control, with performance on a par with larger and more expensive spectrometers.

There are over 45 years of metals industry experience behind the design of the German-made OE series. From measuring nitrogen in steel and iron to phosphorus in aluminium, these analyzers give you the comprehensive metals analysis you need to meet today's tough specifications.

For fast, comprehensive metals quality analysis in a single affordable instrument, the OE series from Hitachi High-Tech delivers everything you need.

Why choose the OE series?

1. Results you can trust

High optical resolution to control tramp and trace elements.

2. Low cost of ownership

Affordable to buy and run, delivering the performance of more expensive analysers.

3. Keeping you running

Engineered for reliability with minimum maintenance and standardization requirements.

4. Supports continuous production

Cloud-based data management, rapid analysis and short start-up times for fast feedback on melt quality.



How we make high-performance affordable

It's all down to the technology

The OE series includes state of the art semiconductor detectors and a new optical concept. This gives high optical resolution, so you don't have to choose between high performance and low cost.

Innovative use of dynamic CMOS detectors and direct coupling of the optics to the spark-stand ensures the best luminosity together with a wavelength range of 174 nm to 670* nm for the OE720 and 119 nm to 670* nm for the OE750.

This covers all elements, from hydrogen to uranium, for complete metals analysis**.

This performance is usually only available with high end quality control instrumentation, yet the OE series keeps costs down through innovation, low argon and power consumption.

The OE series gives you maximum flexibility for element selection, making your operation fit for the future.

Maximum reliability minimum downtime

- The OE series is engineered to keep you running. Maintenance time is minimized; for example you only need to clean the spark stand every few thousand measurements* and standardization intervals are typically weeks, or even months. This could potentially buy you up to two hours of extra productivity.
- Plus, a newly developed electrical spark source gives you better reliability and the best energy excitation pulse per element, making your analysis more precise.
- State-of-the-art CMOS high dynamic detectors

*766 nm on request

**Depending on the application, for further details see our application reports: hhtas.net/stationary-oes



Software innovation to support quality control

The OE series comes with software that makes analysis faster, more accurate and easier to interpret. This offers analysis that helps with process control, and data management functionality allows full traceability of results and helps when it comes to audit time.



SpArcfire

The modern operating software for the OE series with a state-of-the-art user interface is designed to meet the requirements of metallurgical experts but is intuitive enough to be operated by inexperienced users.



GRADE Database

The largest metals database on the market for fast and easy grade identification is pre-installed on the OE series. The Hitachi GRADE Database offers more than 15 million records for over 350,000 materials from over 74 countries and standards. You can update your instrument's grade database with a few clicks – no time-consuming research in norms and grade catalogues.



ExTOPE Connect

ExTOPE Connect is an advanced data management and storage service that allows you to store your results safely, share data instantly and access data in real time from any computer. It includes unlimited free and secure data storage and you can manage a fleet of instruments across several sites from one centralized location.

Technical specifications

OE750
Analyses hydrogen,
nitrogen and oxygen



	OE720	OE750
Dimensions		
Width Height Depth	425 mm 535 mm 760 mm	
Weight	84 kg (185 lbs)	88 kg (194 lbs)
Floor stand version Width Height Depth	425 mm 535 mm 1250 mm	
Floor stand version Weight	162 kg (357 lbs)	166 kg (366 lbs)
Electrics		
Power	100 - 240 V AC, 50 60 Hz	
Consumption max.	430 W	
Operating mode	45 W	
Standby mode	50 W source on	
Optical system		
Rowland circle	Paschen-Runge mounting	
High resolution multi-CMOS	Optimized pixel resolution	
Wavelength range	174 - 670 nm (extendable to 766 nm on request)	119 - 670 nm (extendable to 766 nm on request)
Focal length	400 mm	
Solid state source		
	Computer controlled parameters	
Frequency	80 - 1000 Hz	
Voltage	250 - 500 V High Energy Pre Spark (HEPS)	
Readout system		
External PC workstation	Microsoft® Windows® user interface	
Options	Adapters Spare parts kit Consumables Sample preparation devices Floor stand version	



Reducing operating costs

The OE series punches well above its weight, giving you the performance of an expensive analyzer in an accessible package. We've reduced running costs too. Here's how:

- **Reduced argon consumption**
The OE series uses 10% less argon than its predecessor, thanks to a reduced volume housing the optics and better gas flow around the spark.
- **Low energy consumption**
These analyzers use less than half the energy of previous models, due to an innovative spark generator and smart power management.

Our services

Our global network of service hubs offers a full range of technical support to keep you up and running.

Telephone help-desks

For a fast response to your problem.

Consumables and accessories

From electrodes to spare part kits.

Recertification and maintenance

Ensures your analyser produces the right result year after year.

Training

Understand your analyser and its features.

Extended warranties

Avoid unplanned costs.

Repairs

Fast and efficient turnaround.



Other Products

We've been providing industrial analysis products for the manufacturing industry for 60 years.

Mobile and portable OES

for high performance analysis of alloyed and trace elements; nitrogen analysis in duplex steels.

Handheld XRF

for fast, reliable, non-destructive identification and analysis of alloys.

Handheld LIBS

latest technology for 1-second alloy identification with no X-rays.

More information

To find out more about the OE series visit:
hhtas.net/oe-series

To find out more about our spark spectrometers visit:
www.hitachi-hightech.com/hha

This publication is the copyright of Hitachi High-Tech Analysis Corporation and provides outline information only, which (unless agreed by the company in writing) may not be used, applied or reproduced for any purpose or form part of any order or contract or regarded as the representation relating to the products or services concerned. Hitachi High-Tech Analysis Corporation's policy is one of continued improvement. The company reserves the right to alter, without notice the specification, design or conditions of supply of any product or service. Hitachi High-Tech Analysis Corporation acknowledges all trademarks and registrations.

© Hitachi High-Tech Analysis Corporation, 2025 All rights reserved.

Part number: 10023241
April 2025

HITACHI