

General:	
Pixel size at maximum magnification	0.35 μ m
True low-contrast 3D resolution	2 μ m
Maximum scanning diameter (single/offset)	50 mm / 100 mm
Maximum scanning length	130 mm
Maximum physical object length	200 mm
Size / weight	1260W x 638D x 510H mm / 260 kg
Shielded feedthrough for user's cables/pipes	yes
Power & control signals for optional in-situ stages	yes
Integrated micro-positioning stage	yes
X-RAY SOURCE	
Emitter	Transmission target
Maximum voltage	110 kV
Maximum power	16 W
Smallest spot size	2 μ m
Filter changer, number of positions	7
X-RAY CAMERA	
Camera sensor	Cooled 16 Mp CCD
Camera field of view	72 x 48 mm